



Technique Chart for DVR w/ Sedecal Generator

	1 - 10 cm			11 - 20 cm			21 - 30 cm			31 - 40 cm		
Body Part	KVP Low	KVP Hi	mAs	KVP Low	KVP Hi	mAs	KVP Low	KVP Hi	mAs	KVP Low	KVP Hi	mAs
Abdomen	70	80	5	75	85	8	80	90	10	85	90	12.5
Thorax	75	85	3.2	85	110	3.2	105	120	4	105	120	5
Pelvis	70	75	10	75	85	12.5	80	90	16	80	90	20
Spine	75	80	8	75	85	12.5	80	90	16	80	90	20
Skull	70	75	6.4	75	85	10	75	90	12.5			
Extremities	70	80	3.2	75	85	8						
Exotics	70	80	3.2	75	85	6.4						
Dental	70	80	4	75	85	8						